


















DigitalScope Technical Requirements




DigitalScope™ is a web-based whole slide image (WSI) retrieval and viewing system. It is best viewed with the latest versions of Edge, Chrome, Firefox or Safari. Internet Explorer 11 does not fully implement the requirements for the image viewer, therefore, some features may not work as expected.




| DigitalScope Viewer Tool Functions (<i>Tools in italics are displayed if available</i>) | |
|--|---|
| Icon | Tool |
|  | Navigation - Move the image and change the magnification |
|  | Full Screen - Toggle on and off the full screen function |
|  | <i>Case History</i> - <i>View the case history if available</i> |
|  | <i>Discussion</i> - <i>View the case discussion if available</i> |
|  | <i>Regions of Interest</i> - <i>View the preselected regions of interest if present</i> |
|  | Image Adjustments - Adjust the brightness and contrast |
|  | Overlays - Insert gradient overlays for cell counts |
|  | Tracker -View thumbnail highlighting the current and previously viewed regions |
|  | Measuring Tool - Measure in microns distances between selected points |
|  | Link Tool -Generate links to the image with selected attributes |
|  | <i>Annotations</i> - <i>Allow administrators to annotate images, not visible with routine viewing</i> |
|  | Help - Display additional information |

| Navigation Toolbar Functions | |
|---|--|
| Icon | Tool |
|  | Home - Return to the 0x magnification |

| | |
|---|---|
|  | Panning - Move the image in the chosen direction The image can also be moved using your mouse to click-and-drag |
|  | Magnification Slider - Adjust the magnification The magnification can also be adjusted by scrolling the mouse wheel |
|  | Magnification Buttons - Adjust to a selected magnification |

| Regions of Interest Toolbar Functions | |
|---|--|
| Icon | Tool |
|  | Region of Interest Information - Display information for specific regions of interest |

| Image Adjustments Toolbar Functions | |
|---|--|
| Icon | Tool |
|  | Brightness - Adjust the brightness of the image |
|  | Contrast - Adjust the contrast of the image |
|  | Reset - Reset the brightness and contrast |

| Overlays Toolbar Functions | |
|---|--|
| Icon | Tool |
|  | Miller Disc - Overlays a Miller graticule for cell counts |
|  | Grid - Overlay a grid |
|  | Color Dropdown - Change the color of the overlay |